

















## SYMBOLS GLOSSARY

The following informations are taken from Table 1 - Symbols to convey medical device information of the standard ISO 15223-1:2021 (Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements):

| Reference number and graphic  | Title of symbol               | Description  | Reference number and graphic  | Title of symbol   | Description  |
|---|-------------------------------|--|---|---|--|
| 5.1.1<br>    | Manufacturer                  | Indicates the medical device manufacturer  | 5.1.3<br>    | Date of manufacture   | Indicates the date when the medical device was manufactured  |
| 5.1.4<br>    | Use-by date                   | Indicates the date after which the medical device is not to be used                        | 5.1.5<br>    | Batch code  | Indicates the manufacturer's batch code so that the batch or lot can be identified   |
| 5.1.6<br>  | Catalogue number              | Indicates the manufacturer's catalogue number so that the medical device can be identified | 5.2.4<br> | Sterilized using irradiation                                      | Indicates a medical device that has been sterilized using irradiation  |
| 5.2.6<br>  | Do not resterilize            | Indicates a medical device that is not to be resterilized                                  | 5.2.8<br>  | Do not use if package is damaged and consult instructions for use | Indicates that a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information |
| 5.2.12<br> | Double sterile barrier system | Indicates two sterile barrier systems  | 5.3.7<br>  | Temperature limit   | Indicates the temperature limits to which the medical device can be safely exposed   |

| Reference number and graphic   | Title of symbol | Description  | Reference number and graphic  | Title of symbol                               | Description   |
|--|-----------------|--|---|---|---|
| 5.4.2<br>   | Do not re-use   | Indicates a medical device that is intended for one single use only  | 5.4.3<br>    | Consult instructions for use                  | Indicates the need for the user to consult the instructions for use                                       |
| 5.4.4<br>   | Caution         | Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences | 5.4.8<br>    | Contains biological material of animal origin | Indicates a medical device that contains biological tissue, cells, or their derivatives, of animal origin |
| 5.7.7<br> | Medical device  | Indicates the item is a medical device   | 5.7.10<br> | Unique device identifier                      | Indicates a carrier that contains unique device identifier information                                    |